
```

% A back of envelope calculation to estimate how much signal we will
get
% for the luminometer project, how that will compare to the detector
% noise floor, how much gain we will need, and whether we expect to be
% sensitive enough for the split luciferase assay.

% Paul Lebel, Diane Wiener
% 2020/07/14

clear;

% Physical constants
h = 6.6E-34;
c = 3E8;
lambda_m = 475E-9;
joulesPerPhoton = h*(c/lambda_m);
nAv = 6.02E23;
electronCharge = 1.6E-19;

% Our system, independent parameters
% Transimpedance gain, in Ohms
gain = 1E9;
% Starting point
integrationTime_s = 1;
% Educated guess
collectionEfficiency = 0.3;
% From OPT101 datasheet at 475 nm
ampsPerWatt = 0.23;
% (NEP) From OPT301 plot with 100M gain at 1 Hz, which seems like the
same
% sensor but with an expanded range on the plots. NEP improves with
higher
% gain - we might do better because we plan on using more gain than
100M
noiseEquivalentPower_W = 4E-13;
% From OPT101 datasheet
darkCurrent_A = 2.5E-12;
% ADC bit depth
bitDepth = 24;
% ADC reference voltage
refVoltage = 5;

% Nanoluc is consistently reported as being 150x brighter than
% firefly, which I found a reference that reports 1.6/s with 40% QE.
Units
% are photons/(enzyme second)
enzymePhotonRate = 1.6*0.41*150;

% Reagents details fom Susanna, Jim Wells' lab. Relevant range is 1 pM
to 1
% fM, but adding a decade on either side for comparison.
sampleVolume_L = 35E-6;

```

```

enzymeConcentration_M = logspace(-11,-16, 30);

% Calculated values
Nenzymes = nAv*enzymeConcentration_M*sampleVolume_L;
photonsPerSecond = enzymePhotonRate*Nenzymes;
wattsEmitted = photonsPerSecond*joulesPerPhoton;
wattsCollected = collectionEfficiency*wattsEmitted;
photonsCollected = photonsPerSecond*integrationTime_s;
quantizationLowerLimit = refVoltage/(2^bitDepth);
voltageOut = wattsCollected*ampsPerWatt*gain;
noiseEquivalentPower_V = noiseEquivalentPower_W*ampsPerWatt*gain/
sqrt(integrationTime_s);
photonShotNoise_V =
    sqrt(photonsCollected)*joulesPerPhoton*gain*ampsPerWatt/
integrationTime_s;
darkElectrons = darkCurrent_A*integrationTime_s/electronCharge;
darkElectronsRMS = sqrt(darkElectrons);
darkVoltageRMS = darkElectronsRMS*electronCharge*gain/
integrationTime_s;
sbr_nep = voltageOut./noiseEquivalentPower_V;
snr_phot = voltageOut./photonShotNoise_V;
snr_elec = voltageOut./darkVoltageRMS;
snr_tot = voltageOut./(sqrt(darkVoltageRMS^2 + photonShotNoise_V.^2));

% Plotting
figure('Position',[100,50, 600,700]);
% subplot(2,1,1);
% fill([1E-15,1E-12,1E-12,1E-15],[1E-8, 1E-8, 1E0,1E0],
% [.9,1,.9],'linestyle','none'); hold all;
% set(gca, 'xscale','log');
% set(gca, 'yscale','log');
loglog(enzymeConcentration_M, voltageOut, '-', 'linewidth',2);
hold all;
loglog(enzymeConcentration_M,
    quantizationLowerLimit*ones(numel(enzymeConcentration_M),1), '-', 'linewidth',2);
loglog(enzymeConcentration_M, photonShotNoise_V, '-', 'linewidth',2);
loglog(enzymeConcentration_M,
    noiseEquivalentPower_V*ones(numel(enzymeConcentration_M),1), '-', 'linewidth',2);
loglog(enzymeConcentration_M,
    darkVoltageRMS*ones(numel(enzymeConcentration_M),1), '-', 'linewidth',2);
loglog(enzymeConcentration_M, voltageOut./snr_tot, '-', 'linewidth',2);
hold all;

grid;
ylabel('Voltage out','fontsize',16);
legend('Calculated signal',...
    'ADC lower limit',...
    'Photon Shot Noise',...
    'Noise Equivalent Power (NEP)',...
    'Dark current shot noise',...
    'Total shot noise',...
    'location','northwest');

```


```


title({'Nanoluc detection with OPT101 in high angle collection
      module',[num2str(integrationTime_s) ' s integration']})
xlabel('Enzyme concentration [M]','fontsize',16);

% subplot(2,1,2);
figure;
% fill([1E-15,1E-12,1E-12,1E-15],[1E0, 1E0, 1E8,1E8],
% [.9,1,.9],'linestyle','none'); hold all;
% set(gca, 'xscale','log');
% set(gca, 'yscale','log');
loglog(enzymeConcentration_M, sbr_nep,'linewidth',2); hold all;
loglog(enzymeConcentration_M, snr_phot,'linewidth',2);
loglog(enzymeConcentration_M, snr_elec, 'linewidth',2);
loglog(enzymeConcentration_M, snr_tot, 'linewidth',2);

ylim([1E0,1E8]);
grid;
title('SNR metrics');
ylabel('SNR Metrics','fontsize',16);
xlabel('Enzyme concentration [M]','fontsize',16);
legend('Signal to NEP ratio',...
      'Signal to photon shot noise',...
      'Dark current shot noise',...
      'Total shot noise',...
      'location','northwest');

```





Published with MATLAB® R2019b